

Protex™ PE Non-Surgical Isolation Gown

Model SGL204Z

Are you still trusting your critical PPE to sources on the other side of the globe?

Since its founding in 1916, Shawmut Corporation has advanced materials innovation to improve people's lives. Now, using our rapid product innovation capabilities, we have developed a high-tech manufacturing operation in North America to make Protex™ Advanced Protection isolation gowns.

- Recommended use: non-surgical, minimal fluid exposure procedures
- Easy donning and doffing
- Level 2 compliant to liquid challenges associated with ANSI/AAMI PB70:12

Product Details and Key Features

- AAMI Level 2, Single-Use
- Soft, Comfortable PE Film
- Over-the-Head Donning with Open Back Style
- Integrated Ties at Waist
- No-Sew Welded Exterior Seams
- Thumb Hole Wrist Closure
- PFAS Free
- Color: Blue
- Sizes: L and XL
- Cases of 100 gowns

Reorder Information

| Part # | Size | Color |
|-------------|-------------|-------|
| SGL204Z-BL | Large | Blue |
| SGL204Z-BXL | Extra Large | Blue |

Shawmut Corporation
208 Manley Street
West Bridgewater, MA 02379
shawmutcorporation.com
508-588-3300

© 2021 Shawmut Corporation. All Rights Reserved. INT-Lit# 00107/8 10/21

Protex™

Advanced Protection



To Learn More or Place an Order:

MedicalSales@ShawmutCorporation.com

shawmutcorp.com/safety



shawmut
Advanced Materials™