

Protex™ Shield Non-Surgical Isolation Gown

Model SGL201RC

Are you still trusting your critical PPE to sources on the other side of the globe?

Since its founding in 1916, Shawmut Corporation has advanced materials innovation to improve people's lives. Now, using our rapid product innovation capabilities, we have developed a high-tech domestic manufacturing operation in the USA to make Protex™ Advanced Protection isolation gowns.

- Recommended use: non-surgical, minimal fluid exposure procedures
- Easy donning and doffing
- Level 2 compliant to liquid challenges associated with ANSI/AAMI PB70:12

Product Details and Key Features

- AAMI Level 2, Reusable
- Multi-Layer Synthetic Laminate
- Full-Back Overlap with 360° Coverage
- Sewn In Ties at Waist and Neck
- No-Sew Welded Exterior Seams
- Polyester Knit Cuffs
- PFAS Free
- Gown life: 50 launderings
- Color: Blue
- Sizes: L and XL
- Cases of 48 gowns

Protex™

Advanced Protection



To Learn More or Place an Order:

MedicalSales@ShawmutCorporation.com

shawmutcorp.com/safety

Reorder Information

Part #	Size	Color
SGL201RC-BL	Large	Blue
SGL201RC-BXL	Extra Large	Blue

Shawmut Corporation
208 Manley Street
West Bridgewater, MA 02379
shawmutcorporation.com
508-588-3300

© 2021 Shawmut Corporation. All Rights Reserved. INT-Lit# 00114/4 10/21



shawmut
Advanced Materials™