

# Shawmut Disposable 3-Ply Filtering Face Mask

Model SM3300 | ASTM F2100-19 Level 3 Fluid Resistance

## ASTM F2100-19 Level 3 Fluid Resistant 3-Ply Filtering Face Mask

Shawmut – a world leader in advanced materials – developed the Disposable 3-Ply Filtering Face Mask with safety, comfort, and performance in mind. Designed, engineered, and manufactured in the USA, this disposable 3-ply face mask meets ASTM F2100-19 Level 3 Fluid Resistance and features 3-layer protection, a hydrostatic filtration layer, comfortable earloop construction that secures the mask to your face, and an adjustable nose bridge piece engineered to fit a wide range of face sizes and produce a better fit and all day comfort.

Adhering to our custom design specifications and engineered to exceed typical standards, our custom-formulated Shawmut MB™ meltblown material is an incredibly lightweight but highly efficient and effective filtration media that allows for easy breathing. The mask's inner layer is a soft, non-abrasive non-woven.

Shawmut MB™ is made in our West Bridgewater, MA, facility and many of the other mask components are sourced from nearby to ensure the highest US-spec quality.



### Face Mask Features

- Meets ASTM F2100-19 Level 3 Fluid Resistance
- Made in the USA
- 3-layer protection
- Shawmut MB™ custom-formulated hydrostatic filtration layer
- Lightweight and breathable
- Comfortable elastic earloops
- Adjustable nose piece for better fit
- Size: Adult
- Color: White
- Single-use



## PRODUCT DETAILS

### Shawmut Disposable 3-Ply Filtering Face Mask, Model SM3300

**Style:** Pleated face mask, disposable

**Size:** Adult

**Color:** White

**Number of Layers:** 3-Ply

**Individually Packaged:** No

**Packaging:** 50 masks/box, 24 boxes/case

**UPC:** Box 850026398007, Case 10850026398004

**Country of Origin:** Made in USA

**Standard:** Meets ASTM F2100-19 Level 3 Fluid Resistance

#### Technical Data:



- Product dimensions: 6.75" (L) x 3.75" (W)
- Box dimensions: 7.375" (L) x 3.875" (W) x 3.9375" (H)
- Box weight: 0.40 lbs.
- Case Dimensions: 18.125" (L) x 12.75" (W) x 16.25" (H)
- Case Weight: 10.4 lbs.

#### Storage Conditions and Shelf Life:

Shawmut Disposable 3-Ply Filtering Face Masks have a standard shelf life of 5 years when masks are stored in their original packaging within climatic conditions ranging from -4 °F (-20 °C) to +86 °F (+30 °C).

**Shawmut Corporation**  
208 Manley Street  
West Bridgewater, MA 02379

© 2022 Shawmut Corporation. All Rights Reserved. INT-00198/0 5/22



#### Materials:

- Spunbond non-woven polypropylene
- Meltblown non-woven polypropylene
- Elastic earloop straps
- Flexible nose bridge piece
- All materials are PVC free
- Not made with natural rubber latex

### Contact Shawmut Today to Order Disposable Filtering Face Masks

#### To Learn More or Place an Order:

[MedicalSales@ShawmutCorporation.com](mailto:MedicalSales@ShawmutCorporation.com)  
[shawmutcorporation.com/filtering-face-masks/](http://shawmutcorporation.com/filtering-face-masks/)

**Technical Assistance:** 508-588-3300

#### Reorder Information

Part #	Size	QTY/Box	Boxes/Case
SM3300-24	Adult	50	24